Stolik

Overbed table

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English version on next page

Visually impaired people can find PDF files and other information about our products on our website:

<www.viteacare.com>.

👉 Please contact a specialist dealer if necessary.

A visually impaired user can also ask an accompanying person to read the instructions.

Table of contents

Introduction	14
Product characteristics	15
Intended use.....	15
Contraindications.....	15
Safety of use	15
Product description	16
Assembly.....	17
Technical specifications	17
Reporting technical defects and repairing the device	18
Care and maintenance	18
Cleaning.....	18
Storage.....	18
Disinfection	18
Product life cycle	19
Product disposal.....	19
Symbols	19
Contact	20

Introduction

Thank you for choosing our product. We strongly believe that it will meet your expectations. Please read this document before first use. These instructions are provided to give the necessary information on safe use of the device. Proper adherence to the recommendations will significantly contribute to prolonging the lifetime and aesthetics of the device. Ensure that the “warranty card” document is completed, attached and handed in to the distributor, as this will ensure the best service.

ATTENTION!

When using the device, please always observe the basic precautions according to the below information regarding safety. Read the entire instructions carefully before use. Please retain these instructions for use for future reference.

WARNING!

Improper use can lead to death, serious injury, dangerous conditions or malfunction of the product.

ATTENTION!

Improper use may lead to personal injury and/or damage to the product.

SUGGESTIONS!

Follow the instructions below to keep the product in good working condition.

If the product is used by another person, provide them with a copy of these instructions for use.

Instructions for use refer to the design and specifications of the product at the time of publication. Due to design changes, some of the illustrations and pictures in the instructions for use may not correspond to the product purchased. We reserve the right to make design changes.

The manufacturer disclaims any liability for personal injury or property damage that may result from improper or unsafe use of this product.

Product characteristics

The overbed table has a C-shaped frame, which allows it to be moved freely over the bed. It is equipped with a tabletop on which various items can be placed, and casters with a brake for greater stability. The user can adjust the height and angle of the top to suit their individual needs. The set includes a base frame, a top frame, a top and a mounting butterfly.

Intended use

A functional overbed table that provides the patient with a comfortable workspace directly above the bed, allowing the patient to freely eat, read and take notes during recovery - both in the hospital and at home.

Contraindications

1. Patients without medical indications for the use of a specialised rehabilitation bed.
2. Patients who may not be able to use the bed function safely, even with caregiver supervision.
3. Persons with cardiac problems for whom the Trendelenburg position is medically contraindicated.
4. Patients with spinal injuries who require specialised orthopaedic beds with more precise adjustments and supports.

Medical device in accordance with Regulation (EU) 2017/745 of the European Parliament and Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC (Text with EEA relevance.)

Safety of use

WARNING!

Each "serious incident" related to the product shall be reported to the manufacturer and the competent authority of the Member State in which the user resides.

WARNING!

The product must not be modified in any way; this will void the warranty.

WARNING!

mdh Ltd. is not responsible for improper use of the product, failure to comply with safety rules and for use not in accordance with the intended use.

WARNING!

Please read these instructions for use before using the product and keep them for future reference. Ensure that all information, recommendations and warnings in these instructions for use are fully understood.

- a. The device must not be overturned, thrown or loaded beyond the permissible standard specified in the instructions for the type of device.
- b. The appearance of corrosion points does not affect the safety of use, but reduces the aesthetics. The manufacturer does not allow corroded equipment to be used.
- c. Persons using the device are prohibited to smoke, consume alcohol, take intoxicants and psychotropic drugs.
- d. Do not proceed with the installation of the product before reading these instructions for use.
- e. Any changes to the product made without the knowledge of the manufacturer are the full responsibility of the person making these changes. Changes made without authorization, will void the warranty and invalidate the declaration of conformity with applicable directives.
- f. It is forbidden to use the product in a manner other than that described in these instructions for use.
- g. The manufacturer is not responsible for malfunction or damage to objects and persons caused by improper use of the product.
- h. Users must read these instructions for use carefully, follow the information given and the procedures for using and maintaining this product.
- i. When the patient is not under supervision, the table should not be used and the medical bed should be placed in the lower position with the side rails raised to reduce the risk of falling and injury.
- j. If the table is next to a power cord or other equipment, precautions should be taken to avoid damage to them.
- k. Use and maintain the product in accordance with the instructions contained in this document, and use only spare parts and technical support approved by the manufacturer.
- l. Keep this manual for consultation or training of personnel. In case of sale or use by new users, it should be kept with the product.

Product description


Illustration 1. List of components.

Assembly

Table assembly.

The table has two parts, which have to be connected.

- Place the base with castors on level ground and lock the brakes. Loosen the height adjustment screw (Illustration 26.).
- Install the top element with the table-top by sliding the profile into the profile of the base support, paying attention to the plastic insert, which must be inserted together with the top profile.
- After adjusting the height, tighten the adjustment screw and check that the connection is stable.
- The table is adjustable for the angle of the top and the angle of the base support.
- Adjustment is possible by extending the lever shown in Illustration 27. Once the tilt angle is set, check that the locking pin is in the plate hole and that the connection is stable.


Illustration 2. Table assembly.

Illustration 3. Tilt angle adjustment.

Technical specifications

Base frame dimensions	575 x 510 mm
Worktop dimensions	590 x 390 mm
Bracket angle adjustment	72° / 90°
Tabletop angle adjustment	0° - 90°
Height adjustment	660 - 1030 mm
Maximum table load	45 kg
Wheel diameter	2"
Table weight	8,4 kg

ATTENTION!

Once the table is assembled, check that the height adjustment lock of the table is tightened firmly.

Reporting technical defects and repairing the device

If a defect or damage to the product is discovered during use or maintenance, it must be reported immediately to the distributor or manufacturer.

WARNING!

Under no circumstances should an unqualified person repair or alter the product or its design.

The only person allowed to make any repairs or adjustments is an employee authorised by the manufacturer.

Care and maintenance

If a defect or damage to the product is discovered during use or maintenance, it must be reported immediately to the distributor or manufacturer.

Cleaning

Coatings should be cleaned with a cloth moistened with water. The use of mild cleaning agents for medical devices is permitted. Do not use harsh polishing agents, as this may damage the protective coatings. Regularly remove dirt from moving product parts (hinges, clamping screws, locks).

It is not permitted to use:

- solvents,
- toilet cleaners,
- agents with chlorine content,
- sharp brushes and hard objects,
- corrosive cleaners.

Do not clean the product with pressure or steam units!

Storage

Device should be stored in a dry, ventilated room with a relative humidity of no more than 80%.

Disinfection

If the product is used by many people (e.g. in a nursing home), commercially available disinfectants should be used. The top surface must be cleaned before disinfection. Only tested and approved disinfectants must be used for disinfection.

Information on recognised and approved disinfectants and disinfection methods can be obtained from the national health authority.

When using the product in hospitals and clinics for multiple patients, and the clinic or hospital regulations do not state otherwise, it is recommended to disinfect the hand grips (wipe with spirit) before handing the device to another patient. No contact between the user's uncovered body and the fabric is allowed, unless separate sets of upholstery are used.

**WARNING!**

When using disinfectants, external surfaces may be damaged, which may reduce the long-term functionality of the components. When doing so, observe the manufacturer's instructions for use.

Product life cycle

The expected life of the medical device is 5 years if it is used as intended and if all maintenance and repair activities are performed as required. The life cycle of the medical device also depends on the frequency of use, the operating environment and care of the device. The life cycle of the device can be extended using spare parts. Spare parts are normally available for 5 years after the model is withdrawn from production. The expected life cycle of the medical device is not equivalent to an additional guarantee.

Product disposal

The product must be disposed of in accordance with current local and national regulations. For information on the correct disposal of packaging, metal components, plastic components, neoprene, silicone and polyurethane materials, please contact your local waste disposal company or authorised distributor.

Symbols

	Medical device
	Date of manufacture (printed on the product)
	Manufacturer
	The Distributor
	Conformité Européene (compliance with the European Community)
	Catalog number (to re-order)
	Serial number
	LOT number
	Unique Device Identifier

	See the instructions for use
	Warning
	Possibility of jamming fingers
	Temperature limits
	Protect from humidity
	Indoor use only
	Keep away from sunlight

Contact

Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Poland

Phone No. +48 42 674 83 84, Fax +48 42 636 52 21

www.mdh.pl

www.viteacare.com

We reserve the right to make technical changes due to technical progress.